



"WE MAKE NGS BETTER"

HighPrep Viral DNA & RNA Kit - DX

Efficient Isolation of Nucleic Acids from Viruses

Catalog Nos. HPV-DR20E, HPV-DR96E, HPV-DR96x4E, HPV-DR3840E
Manual Revision 1
WI-72-66


- Viral nucleic acid isolation from whole blood serum, plasma, saliva, and other bodily fluids
- Magnetic bead-based chemistry

Instructions For Use

Contents

Product Description and Process	1
Kit Contents, Shipping and Storage, and Safety Info ...	1
Preparation of Reagents	2
Protocol: 50 µL sample volume (96 well format)	3
Protocol: 200 µL sample volume (96 well format)	5
Protocol: 500 µL sample volume (96 well format)	7
Troubleshooting Guide	9
Ordering and Related Product Information	10

EC REP CEpartner4U
 Esdoornlaan 13, 3951 DB Maarn
 Netherlands
 www.cepartner4u.com

 MagBio Genomics, Inc.
 200 Professional Drive
 Gaithersburg, MD 20879
 USA



For *in vitro* diagnostic procedures.

Information in this document is subject to change without notice.

MAGBIO GENOMICS, INC. DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. TO THE FULLEST EXTENT ALLOWED BY LAW, IN NO EVENT SHALL MAGBIO GENOMICS, INC. BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF, WHETHER OR NOT FORESEEABLE AND WHETHER OR NOT MAGBIO GENOMICS, INC. IS ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

TRADEMARKS

The trademarks mentioned herein are the property of MagBio Genomics, Inc. or their respective owners.

Product Description

The HighPrep Viral DNA & RNA Kit - DX is designed for rapid and reliable isolation of viral nucleic acids from whole blood, serum, plasma, saliva, and other bodily fluids. The kit extracts high quality DNA and RNA that is suitable for direct use in most downstream applications such as amplification and enzymatic reactions. This protocol can be used in manual procedures as well as a guideline for adapting the kit to automated instruments.

This product is intended to be used by qualified and trained laboratory professionals only.

Process

Samples are lysed in a specially formulated lysis buffer. Nucleic acids are bound to the surface of MAG-S1 Particles under proper conditions. Proteins and cellular debris are efficiently washed with few wash steps. Pure RNA and DNA are then eluted in nuclease-free water or low ionic strength buffer. Purified RNA or DNA can be directly used in downstream applications without the need for further purification.

Kit Contents and Storage

HighPrep Viral DNA & RNA Kit - DX Catalog No.*	HPV-DR20E	HPV-DR96E	HPV-DR96x4E	HPV-DR3840E	Storage
Number of Preps	20	96	384	3840	
Viral Lysis Buffer	6 mL	30 mL	120 mL	120 mL x 10	15-25°C
RDW Buffer ¹	6 mL	30 mL	120 mL	120 mL x 10	15-25°C
Nuclease-Free Water	8 mL	35 mL	140 mL	140 mL x 10	15-25°C
Pro K Solution ²	230 µL	1.1 mL	4.4 mL	4.4 mL x 10	2-8°C
NBE ³	180 µL	2 mL	8 mL	8 mL x 10	2-8°C
MAG-S1 Particles ⁴	230 µL	1.1 mL	4.4 mL	4.4 mL x 10	2-8°C

¹Ethanol must be added prior to use. See preparation of reagents section.

*Reagent volumes are based on a 200 µL sample size. Once opened, reagents are usable until the expiration date on the product label. Be sure to close the lid firmly before storing reagents for later use.

Shipping and Storage

- ²Pro K Solution comes in a ready to use solution. Ships at room temperature. Store at 2-8°C.
- ³NBE ships at room temperature. Store at 2-8°C.
- ⁴MAG-S1 Particles ship at room temperature. Store at 2-8°C.

Safety Information

Any consumables, including plates, tubes, etc., used to process samples with infectious or microbial hazards should be disposed of in an appropriate biohazard waste bin. When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). The SDS can be downloaded from the "Product Documents" tab when viewing this product at www.magbiogenomics.com.

Preparation of Reagents

Prepare the following components for each kit before use.

Equipment and Reagents to Be Supplied by the User

- Ethanol (100%)

Catalog No.	Component	Add 100% Ethanol*	Storage
HPV-DR20E	RDW Buffer	4 mL	Room Temp 15-25°C

Catalog No.	Component	Add 100% Ethanol*	Storage
HPV-DR96E	RDW Buffer	20 mL	Room Temp 15-25°C

Catalog No.	Component	Add 100% Ethanol*	Storage
HPV-DR96x4E	RDW Buffer	80 mL	Room Temp 15-25°C

Catalog No.	Component	Add 100% Ethanol*	Storage
HPV-DR3840E	RDW Buffer	80 mL	Room Temp 15-25°C

**Ensure bottle/tube lid is closed tightly when preparing and storing reagents.*

Protocol : 50 μ L sample volume (96 well format)





Equipment and Reagents to Be Supplied by the User

- Ethanol (70%)
- Isopropanol (100%)
- Magnetic separation device for 96-well plate
- 96-well microplates (U or V bottom)

Things to do Before Starting

- Prepare all reagents according to the instructions on page 2

Protocol

-  Bring the **MAG-S1 Particles** to room temperature for at least 30 min before use.
1. Add 50 μ L of bodily fluid or cell culture lysate to each sample well.
Note: If sample volume is less than 50 μ L, bring volume up to 50 μ L with **Nuclease-Free Water**.
 2. To 50 μ L of sample add 60 μ L of **Viral Lysis Buffer** and 5 μ L of **Pro K Solution**. Mix very well.
 3. Incubate at 56-60°C for 10 min. A thermoshaker may be used during incubation. If a thermoshaker is not used, shake the samples once or twice during the incubation period.
 4. Let the samples cool to room temperature and add 4 μ L of **NBE**, 70 μ L of 100% Isopropanol, and 5 μ L of **MAG-S1 Particles**. Pipette mix 15 times.
 *Shake well to resuspend the **MAG-S1 Particles** before use.*
 5. Let the samples sit at room temperature for 10 min.
 6. Place the sample plate on the magnetic separation device for 5 min to magnetize the **MAG-S1 Particles** or until the magnetic beads clear from the solution.
 7. With the plate on the magnetic separation device, remove and discard the supernatant by pipetting.
 *Do not disturb the attracted beads while aspirating the supernatant.*
 8. Remove the plate from the magnetic separation device.
 9. Add 200 μ L of **RDW Buffer** to each sample and pipette mix 15 times to resuspend the **MAG-S1 Particles**.
 10. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
 11. With the plate on the magnetic separation device, remove and discard the supernatant.
 *Do not disturb the attracted beads while aspirating the supernatant.*

12. Remove the plate from the magnetic separation device.
13. Add 250 μ L of 70% Ethanol to the sample and pipette mix 15 times to resuspend the **MAG-S1 Particles**.
*⚠ Complete resuspension of the **MAG-S1 Particles** is crucial for obtaining quality DNA or RNA.*
14. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
15. With the plate on the magnetic separation device, remove and discard the supernatant.
⚠ Do not disturb the attracted beads while aspirating the supernatant.
16. Repeat steps 12-15 for a second wash.
17. Dry the beads by incubating for 7 min at room temperature with the plate still on the magnetic separation device.
⚠ It is critical to completely remove any residual liquid from each well.
18. Remove the plate from the magnetic separation device. Add 20-50 μ L of **Nuclease-Free Water** (pre-warmed to 60-65°C) to each well and pipette mix 25 times to completely resuspend the **MAG-S1 Particles**.
*⚠ Complete resuspension of the **MAG-S1 Particles** is crucial for obtaining purity.*
19. Incubate at room temperature for 10 min.
⚠ Incubation at 65°C for 5 min may increase DNA yield, but may affect the quality of RNA.
20. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
21. Transfer the eluate (cleared supernatant containing the DNA or RNA) to a new microplate for storage. Store DNA at -20°C and RNA at -80°C.

Protocol : 200 μ L sample volume (96 well format)

Equipment and Reagents to Be Supplied by the User

- Ethanol (70%)
- Isopropanol (100%)
- Magnetic separation device for 96-well plate
- 96-well microplates (U or V bottom)

Things to do Before Starting


- Prepare all reagents according to the instructions on page 2

Protocol

 Bring the **MAG-S1 Particles** to room temperature for at least 30 min before use.

1. Add 200 μ L of bodily fluid or cell culture lysate to each sample well.
Note: If sample volume is less than 200 μ L, bring volume up to 200 μ L with **Nuclease-Free Water**.
2. To 200 μ L of sample add 240 μ L of **Viral Lysis Buffer** and 10 μ L of **Pro K Solution**. Mix very well.
3. Incubate at 56-60°C for 10 min. A thermoshaker may be used during incubation. If a thermoshaker is not used, shake the samples once or twice during the incubation period.
4. Let the samples cool to room temperature and add 8 μ L of **NBE**, 280 μ L of 100% Isopropanol, and 10 μ L of **MAG-S1 Particles**. Pipette mix 15 times.

 *Shake well to resuspend the **MAG-S1 Particles** before use.*

5. Let the samples sit at room temperature for 10 min.
6. Place the sample plate on the magnetic separation device for 5 min to magnetize the **MAG-S1 Particles** or until the magnetic beads clear from the solution.
7. With the plate on the magnetic separation device, remove and discard the supernatant by pipetting.
 *Do not disturb the attracted beads while aspirating the supernatant.*
8. Remove the plate from the magnetic separation device.
9. Add 400 μ L of **RDW Buffer** to each sample and pipette mix 15 times to resuspend the **MAG-S1 Particles**.
10. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.

11. With the plate on the magnetic separation device, remove and discard the supernatant.

 *Do not disturb the attracted beads while aspirating the supernatant.*

12. Remove the plate from the magnetic separation device.
13. Add 500 μL of 70% Ethanol to the sample and pipette mix 15 times to resuspend the **MAG-S1 Particles**.
*⚠ Complete resuspension of the **MAG-S1 Particles** is crucial for obtaining quality DNA or RNA.*
14. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
15. With the plate on the magnetic separation device, remove and discard the supernatant.
⚠ Do not disturb the attracted beads while aspirating the supernatant.
16. Repeat steps 12-15 for a second wash.
17. Dry the beads by incubating for 10 min at room temperature with the plate still on the magnetic separation device.
⚠ It is critical to completely remove any residual liquid from each well.
18. Remove the plate from the magnetic separation device. Add 50-100 μL of **Nuclease-Free Water** (pre-warmed to 60-65°C) to each well and pipette mix 25 times to completely resuspend the **MAG-S1 Particles**.
*⚠ Complete resuspension of the **MAG-S1 Particles** is crucial for obtaining purity.*
19. Incubate at room temperature for 10 min.
⚠ Incubation at 65°C for 5 min may increase DNA yield, but may affect the quality of RNA.
20. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
21. Transfer the eluate (cleared supernatant containing the DNA or RNA) to a new microplate for storage. Store DNA at -20°C and RNA at -80°C.

Protocol : 500 μ L sample volume (96 well format)





Equipment and Reagents to Be Supplied by the User

- Ethanol (70%)
- Isopropanol (100%)
- Magnetic separation device for 96-well plate
- 96-well deep plate

Things to do Before Starting

- Prepare all reagents according to the instructions on page 2

Protocol

-  Bring the **MAG-S1 Particles** to room temperature for at least 30 min before use.
- 1. Add 500 μ L of bodily fluid or cell culture lysate to each sample well.
Note: If sample volume is less than 500 μ L, bring volume up to 500 μ L with **Nuclease-Free Water**.
- 2. To 500 μ L of sample add 500 μ L of **Viral Lysis Buffer** and 20 μ L of **Pro K Solution**. Mix very well.
- 3. Incubate at 56-60°C for 10 min. A thermoshaker may be used during incubation. If a thermoshaker is not used, shake the samples once or twice during the incubation period.
- 4. Let the samples cool to room temperature and add 16 μ L of **NBE**, 500 μ L of 100% Isopropanol, and 10 μ L of **MAG-S1 Particles**. Pipette mix 15 times.
 *Shake well to resuspend the **MAG-S1 Particles** before use.*
- 5. Let the samples sit at room temperature for 10 min.
- 6. Place the sample plate on the magnetic separation device for 10 min to magnetize the **MAG-S1 Particles** or until the magnetic beads clear from the solution.
- 7. With the plate on the magnetic separation device, remove and discard the supernatant by pipetting.
 *Do not disturb the attracted beads while aspirating the supernatant.*
- 8. Remove the plate from the magnetic separation device.
- 9. Add 500 μ L of **RDW Buffer** to each sample and pipette mix 15 times to resuspend the **MAG-S1 Particles**.
- 10. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
- 11. With the plate on the magnetic separation device, remove and discard the supernatant.
 *Do not disturb the attracted beads while aspirating the supernatant.*

12. Remove the plate from the magnetic separation device.
13. Add 500 μL of 70% Ethanol to the sample and pipette mix 15 times to resuspend the **MAG-S1 Particles**.
*⚠ Complete resuspension of the **MAG-S1 Particles** is crucial for obtaining quality DNA or RNA.*
14. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
15. With the plate on the magnetic separation device, remove and discard the supernatant.
⚠ Do not disturb the attracted beads while aspirating the supernatant.
16. Repeat steps 12-15 for a second wash.
17. Dry the beads by incubating for 10 min at room temperature with the plate still on the magnetic separation device.
⚠ It is critical to completely remove any residual liquid from each well.
18. Remove the plate from the magnetic separation device. Add 50-100 μL of **Nuclease-Free Water** (pre-warmed to 60-65°C) to each well and pipette mix 25 times to completely resuspend the **MAG-S1 Particles**.
*⚠ Complete resuspension of the **MAG-S1 Particles** is crucial for obtaining purity.*
19. Incubate at room temperature for 10 min.
⚠ Incubation at 65°C for 5 min may increase DNA yield, but may affect the quality of RNA.
20. Place the sample plate back on the magnetic separation device and wait for 5 min or until the magnetic beads clear from the solution.
21. Transfer the eluate (cleared supernatant containing the DNA or RNA) to a new microplate for storage. Store DNA at -20°C and RNA at -80°C.

Troubleshooting Guide

Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact technical support via:

Phone: US/Canada, +1 301-302-0144. Europe, +49 7250 33 13 403

Email: US/Canada, support@magbiogenomics.com. Europe, Info.europe@magbiogenomics.com

Symptoms	Possible Causes	Comments
Low DNA or RNA yield	Incomplete resuspension of MAG-S1 Particles	Resuspend MAG-S1 Particles by vortexing vigorously before use
	Loss of MAG-S1 Particles during operation	Avoid disturbing the MAG-S1 Particles during aspiration of supernatant
	Ethanol is not added into RDW Buffer	Add absolute 100% Ethanol to RDW Buffer (see page 2 for instructions)
	Inefficient cell lysis	Double the volume of Pro K Solution and incubate longer
Low sample purity	Inefficient sample wash	Perform an additional wash step with RDW Buffer
MAG-S1 Particles do not completely clear from the solution	Magnetizing time too short	Increase collection time on the magnet
Problems in downstream applications	Insufficient DNA/RNA in starting material	Use more starting material.
	Ethanol carry-over	Dry the MAG-S1 Particles completely before elution
Carryover of MAG-S1 Particles	The eluate has particles and is not fully clear	Increase magnetization time. If small amount of carryover, place eluted sample on a magnetic separation device and perform an additional 5 min magnetization

Ordering Information

HighPrep Viral DNA & RNA Kit - DX

Catalog No.	Product	Description	Preps
HPV-DR96E	HighPrep Viral DNA & RNA Kit - DX (96 Preps)	Magnetic bead-based kit for rapid and reliable isolation of viral nucleic acids from whole blood, serum, plasma, saliva, and other bodily fluids.	96
HPV-DR96x4E	HighPrep Viral DNA & RNA Kit - DX (384 Preps)	Magnetic bead-based kit for rapid and reliable isolation of viral nucleic acids from whole blood, serum, plasma, saliva, and other bodily fluids.	384
HPV-DR3840E	HighPrep Viral DNA & RNA Kit - DX (3840 Preps)	Magnetic bead-based kit for rapid and reliable isolation of viral nucleic acids from whole blood, serum, plasma, saliva, and other bodily fluids.	3840

Related Products

HighPrep PCR - DX

Catalog No.	Product
AC-60005E	HighPrep PCR - DX (5 mL)
AC-60050E	HighPrep PCR - DX (50 mL)
AC-60250E	HighPrep PCR - DX (250 mL)
AC-60500E	HighPrep PCR - DX (500 mL)

HighPrep RNA Elite - DX

Catalog No.	Product
RC-90005E	HighPrep RNA Elite - DX (5 mL)
RC-90050E	HighPrep RNA Elite - DX (50 mL)
RC-90250E	HighPrep RNA Elite - DX (250 mL)
RC-90500E	HighPrep RNA Elite - DX (500 mL)

Magnetic Separation Devices

Catalog No.	Product
MYMAG-96	Handheld Magnetic Separation Device (96 well microplate format)
MYMAG-96X	Magnetic Separation Device (96 well ring format)
MBMS-12	MagStrip Magnet Stand (1.5 mL x 12)
MBMS-31550	15 mL and 50 mL Magnetic Stand Combo (3 x 15 mL and 3 x 50 mL)



www.magbiogenomics.com